ENGLISH

A. DEVICE DESCRIPTION

The WasherCap™ Mini Implantable device consists of a cap and screw for suture/tape fixation to bone. The WasherCap™ Mini Fixation System consists of the WasherCap™ Mini Implantable device, an inserter, a screwdriver, a drilling guide, a drill bit and a suture threader.

The WasherCap™ Mini Implantable device is pre-loaded on the inserter with the suture threader. The drilling guide and the drill bit should be used to perform a drilling guided specifically to the angulation necessary to create the bone socket where the screw of the implant should be located. The suture/tape is passed through the cap with the suture threader and the screw is threaded with the screwdriver, setting the tension of the suture/tape.

The WasherCap™ Mini Fixation System is available in different sizes: 3.5mm and 4.5mm.

B. INDICATIONS FOR USE

WasherCap™ Mini Fixation System is intended for fixation of suture/tape (soft tissue) to bone in the shoulder, knee and hand/wrist, in skeletally mature pediatric and adult patients for the following procedures:

- **Knee:**
 - **Meniscal root repair:** using non-absorbable UHMWPE USP 2-0 sutures or non-absorbable UHMWPE 1.4-2.2mm tapes.
 - **Meniscal transplant:** using non-absorbable UHMWPE USP 2-0 sutures or non-absorbable UHMWPE 1.4-2.2mm tapes.
 - **ACL/PCL Repair:** using non-absorbable UHMWPE USP 2-0 sutures or non-absorbable UHMWPE 1.4-2.2mm tapes.
 - **ACL Reconstruction: only with WasherCap™ Mini Fixation System 4.5mm size and using 2 non-absorbable UHMWPE 1.4mm tapes or 1 non-absorbable UHMWPE 2.2mm tape.**
- **Shoulder:**
 - **Rotator Cuff Repair:** using non-absorbable UHMWPE USP 2-0 sutures or non-absorbable UHMWPE 1.4-2.2mm tapes.
- **Hand/Wrist:**
 - **TFCC: only with WasherCap™ Mini Fixation System 3.5mm size and using non-absorbable UHMWPE USP 2-0 sutures or non-absorbable UHMWPE 1.4-2.2mm tapes.**

IMPORTANT NOTE: See the WasherCap™ Mini Fixation System Surgical Technique for the correct method of tape and/or suture fixation.

NOTE: Maximum capacity of sutures or tapes per size (*):

Size	2-0 suture	1.4mm tape	2mm tape	2.2mm tape
WasherCap™ Mini 3.5mm	7 sutures (14 suture ends)	2 tapes (4 tape ends)	1 tape (2 tape ends)	1 tape (2 tape ends)
WasherCap™ Mini 4.5mm	9 sutures (18 suture ends)	3 tapes (6 tape ends)	2 tapes (4 tape ends)	1 tape (2 tape ends)

(*): Ensure that your configuration fits through the device correctly before performing surgery.

C. CONTRAINDICATIONS

- The use of WasherCap™ Mini Fixation System is contraindicated in the following cases:
- Limited blood supply and previous infections which may retard healing.
 - Foreign body sensitivity and known hypersensitivity to the implant materials. Where foreign body and material sensitivity is suspected, appropriate tests should be carried out and sensitivity ruled out prior implantation.
 - Any active infections, either local or systemic.
 - Conditions which may limit the patient's ability to follow a post-operative regimen or willingness to restrict activities during the healing period.
 - Insufficient quantity or quality of bone.
 - The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing the surgery on patients who are skeletally immature.
 - Do not use the WasherCap™ Mini Fixation System for other surgeries than the indicated in these Instructions For Use.
 - Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.
 - Do not use sutures and tapes manufactured with materials not indicated in these instructions for use. Only FDA certified non-absorbable UHMWPE or non-absorbable UHMWPE + PET tracer's sutures and tapes can be used.
 - Any concomitant disease that may put at risk the stability of the implant upon doctor's judgment.

D. ADVERSE EVENTS

- Infections, both deep and/or superficial.
- Allergies and other reactions to the device materials including foreign bodies reactions.
- Implant failure leads to removal and follow-up surgeries.

E. WARNINGS

- Caution: Federal (US) law restricts this device to sale by or on the order of a physician.
- WasherCap™ Mini Fixation System should only be implanted by trained surgeons who are familiar with the clinical use, surgical techniques, proper selection and placement of the WasherCap™ Mini Implant device and its associated accessories for achieving a good surgical result and minimizing the risks to patients.
- Size selection of the implant should be made with care taking into consideration the indications of use included in these instructions for use.
- WasherCap™ Mini Fixation System (implantable device and its associated single use surgical accessories) is provided STERILE and is a SINGLE PATIENT USE device. Do NOT RE-STERILIZE, do not re-use or re-implant this device under any circumstances. Re-use of the WasherCap™ Mini Fixation System could cause risks of malfunction, breaking which may compromise device performance as intended and could cause harm to the patient and/or user. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- Preoperative and surgical preparation should be evaluated based on knowledge of surgical techniques and the device (including selection of the appropriate size according to its intended use, its implantation and the use of associated surgical accessories) to avoid risks to the patient associated with the duration of the operation and to ensure a satisfactory surgical result and use of the device.
- During the postoperative period until healing is complete the fixation provided by this device should be protected and considered temporary as it cannot support weight bearing or other unsupported stresses. The postoperative regimen prescribed by the physician should be strictly followed. Failure to comply with the postoperative regimen could lead to device failure and affect the expected results.
- After surgery, wounds or lesions on the site of intervention could appear and/or pain, discomfort or abnormal feeling due to the presence of the device.
- Any decision to remove the device should be evaluated based on the risk posed to the patient by a second surgical intervention. If the device is removed, adequate postoperative management should be followed.
- Ensure that the associated surgical accessories included in the WasherCap™ Mini Fixation System are used to use the implant. This device.
- If any component is damaged, a new package of the WasherCap™ Mini Fixation System of the same size must be opened to obtain a new component. If the size of any of the components does not match the rest of the components, the integrity of the other components may be compromised. In the case of the implant, a size that does not match the size of the bone tunnel can jeopardize the fixation.
- The screwdriver included in the WasherCap™ Mini Fixation System is designed to be a torque limiter and break (without any risk) if it is over-tightened in order to maintain the surgical security. In the event that the screwdriver breaks, and the user wishes to adjust the suture and tape fixation, a new package of the WasherCap™ Mini Fixation System of the same size must be opened to obtain a new screwdriver.
- Biohazard waste, such as explanted devices, single use surgical accessories, and contaminated surgical equipment, should be safely disposed of in accordance with the institution's policy.
- Serious incidents involving the device should be reported to ABANZA or a representative in the country, and to the health authorities corresponding to the place where the incident happened.
- In the case the implant requires removal and a follow-up surgery is scheduled, a device of the same size or one size larger than the one initially implanted should be used. In any case, please contact your ABANZA representative or your local distributor.

F. MRI (MAGNETIC RESONANCE IMAGING) SAFETY INFORMATION

The WasherCap™ Mini Fixation System has not been evaluated for safety in the MR environment. It has not been tested for heating or unwanted movement in the MR environment. The safety of WasherCap™ Mini Fixation System in the MR environment is unknown. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

G. PRECAUTIONS

- When determining the appropriate size of the WasherCap™ Mini Fixation System, surgeons should use professional judgment based on the patient's medical history, the surgical technique to be employed and the indications for use.
- Prior to performing any surgery with the WasherCap™ Mini Fixation System, surgeons are advised to review the device-specific surgical techniques. ABANZA offers a detailed description of surgical techniques in a variety of formats (print, video, electronic). More detailed information on surgical techniques and demonstrations can be found on the ABANZA website. For an on-site demonstration, please contact your ABANZA representative.
- Measurement of the tunnel diameter and the proper implant selection are critical to maintain functionality of the implant construct.
- The cap must be properly oriented during insertion in order to function appropriately.
- Inserting the cap with an excessive angle than the recommended one may damage the device or suture construct. Inserting the WasherCap™ Mini Implantable device in an incorrect position in the bone tunnel axis may result in implant failure or excessive protrusion from the cortical bone.
- Multiple insertions of the device may damage the device or suture construction.
- During the screw insertion phase, it is important for the user to use only rotary movements of the screwdriver and not to combine them with an axial insertion force. An axial insertion force on the screw may cause implant failure. See Surgical Techniques for more information (available on the ABANZA website). ABANZA recommends that the screw should not be fully unscrewed to re-tension the sutures and tapes. Completely unscrewing the screw from its cap may cause complications during reassembly and fixation. If the screw is completely removed from the cap, ABANZA recommends the use of a new device of the same size to ensure the safety of the implant.
- Do not use the implant device or associated surgical accessories if it gets contaminated (in contact with non-sterile surfaces) or suspected to be contaminated with other elements prior to or during the surgery.
- Some instruments included in the WasherCap™ Mini Fixation System, due to their intended use, have sharp edges that can cause damage on sutures, tapes, patient or user.

H. PACKAGING AND LABELLING

Carefully inspect the device before use. Do not use the device if you suspect that the package has been opened and/ or that any component of the device is damaged. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature, or has been altered in any way. Please contact ABANZA Customer Service to report any package damage or alterations.

Contact customer service at "customerservice@abanzamed.com".

I. STERILIZATION

WasherCap™ Mini Fixation System (implantable device and its associated single use surgical accessories) is provided STERILE and ready to use. Refer to the package label for the sterilization method. WasherCap™ Mini Fixation System must not be re-sterilized under any circumstance as it is provided in a terminally sterilized configuration.

J. MATERIAL SPECIFICATIONS

- **WasherCap™ Mini implantable device** – the implantable part of the WasherCap™ Mini Fixation System is composed of:
 - Cap: PEEK – Polyetheretherketone and titanium alloy.
 - Screw: Titanium alloy.
- **WasherCap™ Mini associated surgical accessories:**
 - **Inserter – WasherCap™ Mini:** this component is for single use and is not intended to be implanted.
 - Metal part – Stainless Steel. This part of the inserter is for single use and intended for transient use.
 - Handle part – Polyamide (PA). This part of the inserter is for single use and intended for non-invasive use.
 - **Screwdriver – WasherCap™ Mini:** this component is for single use and is not intended to be implanted.
 - Metal part – Stainless Steel. This part of the screwdriver is for single use and intended for transient use.
 - Handle part – Polyamide (PA). This part of the screwdriver is for single use and intended for non-invasive use.
 - **Suture threader – WasherCap™ Mini:** Nitinol and Stainless Steel. This component is for single use and intended for non-invasive use (not intended to be implanted).
 - **Drilling guide – WasherCap™ Mini:** Polyamide (PA). This component is for single use and intended for transient use (not intended to be implanted).
 - **Drill bit – WasherCap™ Mini:** Stainless Steel. This component is for single use and intended for transient use (not intended to be implanted).

K. STORAGE CONDITIONS

- WasherCap™ Mini Fixation System must be stored within its original unopened packaging.
- WasherCap™ Mini Fixation System should be stored in a dry place (away from moisture) away from sunlight and must not be used beyond the expiration date indicated on the package.

L. DISCLAIMER

Please contact ABANZA's customer service for a copy of the Instructions For Use (IFU). All implant devices and/or the associated single use surgical instrumentation which may be damaged, scratched, showing defects, corrosion, cracked joints or dislocation, should not be used and be managed according to hospital protocol rules. Any misuse, modification or substitution of the WasherCap™ Mini Fixation System It is strictly forbidden. If this recommendation is not followed by the surgical team, ABANZA rejects any responsibility derived from any subsequent consequence.

This document is not subject to control after distribution, and thus, it is the responsibility of the reader to confirm that this is the most current revision or to contact ABANZA's customer service to request the latest version.

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