WASHERCAP™
FIXATION SYSTEM

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ENGLISH

A. DEVICE DESCRIPTION

WasherCap™ Fixation System is a device which provides soft tissue graft fastening, available in a range of sizes. WasherCap™ Fixation System consists of three parts:

- CAP**: main component of the System where the graft compression is achieved. It is made of biocompatible polyetheretherketone (PEEK).
- WASHER**: this part is provided together with the Screw and it is made of biocompatible PEEK. It should never be separated from the Screw.
- SCREW**: made of surgical biocompatible titanium alloy. It should never be separated from the Washer. Sizes: the different sizes are based on the thickness of the graft used during the surgery procedure. The following sizes are available: 7mm, 8mm, 9mm, 10mm.

WasherCap™ Fixation System should be implanted by a trained trauma surgeon following the surgical procedure recommended in the given Instructions for Use (IFU).

As it is stated in section "E. PRECAUTIONS", please carefully inspect the device package prior implantation to ensure it has not been opened and/or any component of the device is damaged. In case of the product has been compromised or damaged, please contact your ABANZA representative or your local distributor.

B. INDICATIONS FOR USE

WasherCap™ Fixation System is intended for fixation of soft tissue grafts, including tendons and ligaments, during surgical procedures such as in Anterior Cruciate Ligament (ACL) reconstruction of the knee.

C. CONTRAINDICATIONS

The use of WasherCap™ Fixation System is contraindicated in the following cases:

- Fixation using Bone-tendon-bone graft.
- Known hypersensitivity to the implant materials. Where material sensitivity is suspected, appropriate tests should be carried out and sensitivity ruled out prior implantation.
- The use of this device may not be suitable for patients with **insufficient or immature bone**. The physician should carefully assess **bone quality** before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Conditions which may limit the patient's ability to follow a post-operative regimen or willingness to restrict activities during the healing period.
- Active infections, either local or systemic.
- Limited blood supply and previous infections which may tend to retard healing.
- Any concomitant disease that may put at risk the stability of the implant upon doctor's judgement (for example, severe osteoporosis).

D. WARNINGS

- WasherCap™ Fixation System is provided **STERILE** and is a single patient use device. Do NOT RE-STERILIZE, do not re-use or re-implanted this device under any circumstances. Re-use WasherCap™ Fixation System could cause risks of malfunction, breaking or infection. Do not use after the expiration date.
- WasherCap™ Fixation System should only be implanted by trauma surgeons who are familiar with the clinical use, surgical technique, proper selection and placement of the device and use of the WasherCap™ associated instrumentation.
- Federal (US) law restricts this device to sale by or on the order of a physician.
- Instruments are provided non-sterile and must be cleaned and sterilized prior to use following the "Reprocessing Instructions" (CI-001).
- Do not use the product for other purposes or surgical procedures other than those indicated in this Instructions for Use.
- In case of removal, the device may result in a potential biohazard and should be handled and disposed in accordance with accepted medical practices as well as applicable local and national requirements.

E. PRECAUTIONS

- Carefully inspect the device before use. Do not use the device if you suspect that the package has been opened and/or that any component of the device is damaged.
- Ensure the device is totally implanted to avoid possible discomfort after the surgery.
- Any decision to remove the implanted device should take into consideration the potential risk to the patient of a second surgical intervention. Implant removal should be followed by adequate postoperative management.
- Exercise caution when considering the use of MRI. The WasherCap™ Fixation System has not been evaluated for safety and compatibility in the MR environment and may cause interference as it has not been tested. The WasherCap™ Fixation System has not been tested for heating or migration in the MR environment.

Note: when choosing the appropriate size of the device, the surgeon will need to consider its intended use, the surgical technique to be used, and the patient's medical history to make a decision based on their professional judgment.

F. POTENTIAL ADVERSE EVENTS

- Allergies and other reactions to the device materials including foreign bodies.
- Infections, both deep and/or superficial.
- Implant failure leading removal and follow-up surgeries.
- Loss of fixation or implant breakage attributable to delayed or non-union.
- Disassembly, bending, fracture, loosening or migration of the implant.
- Decrease of bone density due to stress shielding.
- Fracture of bony structures.
- Vascular damage resulting in catastrophic bleeding and/or respiratory complications, including pulmonary effusion and pneumothorax and/or nerve damage, including loss of neurologic function due to surgical trauma or presence of the device.
- Death (due to clinical intervention).

These adverse effects include adverse effects that are important considerations for internal fixation devices.


G. STERILIZATION

- WasherCap™ Fixation System is provided **STERILE** and ready to use. Refer to the package label for the sterilization method. The implant must not be re-sterilized under any circumstance.
- The associated **Surgical Instrument Set** is not provided sterile. All surgical instruments must be cleaned and sterilized before use as detailed in the "Reprocessing Instructions" (CI-001). Instruments may be re-used after cleaning and sterilization.
- Recessed and hidden areas within instruments should be inspected to ensure that entrapped or other residual materials are completely removed.
- The sterilization parameters described in the "Reprocessing Instructions" (CI-001) are only valid for surgical instrument sets that are adequately cleaned.
- It is recommended that the associated **Surgical Instrument** be sterilized only by moist heat autoclaving procedures regularly used in the hospital, as detailed in the "Reprocessing Instructions" (CI-001).
- The document "Reprocessing Instructions" (CI-001) indicates which instruments must be disassembled before cleaning and reassembled before sterilization.

Note: the autoclave must be properly installed, maintained, and calibrated. The autoclave manufacturer's operating instructions and recommended guidelines for maximum sterilization load should be followed. Your institution is responsible for the normal sterility assurance validation.

This document is not subject to control after distribution, and thus, it is the responsibility of the reader to confirm that this is the most current revision or to contact ABANZA TECNOMED's customer service to request the latest version.

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