ENGLISH

These instructions were developed using the guidance given in the following standards:

- ISO 17664: Processing of health care products – Information to be provided by the medical manufacturer for the processing of medical devices.
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
- ANSI/AAMI ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in health care Facilities.
- ISO 16061: Instrumentation for use in association with non-active surgical implants – general requirements.
- Guidance for Industry and Food and Drug Administration Staff: "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling". June 9, 2017.

A. INSTRUMENT SET INFORMATION AND DESCRIPTION

Intended use:

The WasherCap™ Surgical Instrument Set is a reusable surgical instrumentation intended to be used in the context of a planned Anterior Cruciate Ligament (ACL) reconstruction surgery, carried out with the WasherCap™ Fixation System with the aim of ensuring that the surgical procedure and the use of the WasherCap™ are performed correctly. This set is suitable for all available models (19 & 21) and sizes of the WasherCap™ Fixation System. The WasherCap™ Surgical Instrument Set is composed of reusable surgical instruments, made of surgical stainless steel according to ASTM F-899 and ASTM F138.

The Set ABWCIN00 is composed of:

- 1 Sterilization tray.
- 4 Countersink Positioners for the complete range of the WasherCap™ device (7mm, 8mm, 9mm and 10mm).
- 1 Positioner guide.
- 2 Kirschner wires.
- 4 Countersinks (consisting of 2 pieces: 1 Main body and 1 Stop element) for the complete range of the WasherCap™ device (7mm, 8mm, 9mm and 10mm).
- 1 Handle with AO connection system.
- 2 Inserters for the complete range of the WasherCap™ device (7mm & 8mm and 9mm & 10mm).
- 2 Inserter Coupling Screws for the complete range of the WasherCap™ device (7mm & 8mm and 9mm & 10mm).
- 1 WasherCap™ Screwdriver.
- 1 Countersink tool.

The Set ABWCIN01 is composed of:

- 1 Sterilization tray.
- 4 Countersink Positioners for the complete range of the WasherCap™ device (8mm, 9mm, 10mm and 11mm).
- 1 Positioner guide.
- 2 Kirschner wires.
- 4 Countersinks (consisting of 2 pieces: 1 Main body and 1 Stop element) for the complete range of the WasherCap™ device (8mm, 9mm, 10mm and 11mm).
- 1 Handle with AO connection system.
- 3 Inserters for the complete range of the WasherCap™ device (8mm, 9mm & 10mm and 11mm).
- 3 Inserter Coupling Screws for the complete range of the WasherCap™ device (8mm, 9mm & 10mm and 11mm).
- 2 Awls for the complete range of the WasherCap™ device (8mm and 9mm & 10mm & 11mm).
- 1 WasherCap™ Screwdriver.
- 1 Countersink tool.

Those instruments available in various sizes can be identified by the size laser marking and/or by the colored marking on each instrument. Each color represents a size: light blue (7mm), yellow (8mm), red (9mm), green (10mm) and dark blue (11mm). ABANZA provides detailed surgical techniques and detailed reprocessing instructions in print, video and electronic formats. ABANZA also provides demonstrations and information on detailed surgical techniques and detailed reprocessing instructions (www.abanzamed.com).

B. LIMITATIONS ON REPROCESSING

Repeated reprocessing has minimal effect on the components of the WasherCap™ Surgical Instrument Set when performed according to ABANZA's "Reprocessing Instructions" (CI-001). End of useful life is normally determined by wear and damage due to use. In addition, it is recommended that, prior to re-using the components of the Surgical Instrument Set (after cleaning process), a visual inspection is carried out to verify that no signs of visible dirt, corrosion and/or other type of damage are present. Please contact your local distributor or ABANZA representative for more information on the reuse test and if you detect wear and/or damage on the components of the Surgical Instrument Set for replacement.

C. REPROCESSING VALIDATION

When properly performed cleaning, disinfection (if applicable) and sterilization, do not compromise the use and mechanical performance of this Surgical Instrument Set. The instruments are used with or on patients who may be harboring either recognized and unrecognized infections. To prevent the spread of infection, all reusable components of the Surgical Instrument Set must be thoroughly cleaned, disinfected (if applicable), and sterilized after use in each patient, following the specifications of the "Reprocessing Instructions" (CI-001) and these Instructions for Use. The recommended reprocessing methods in these IFU have been validated in compliance with federal and international guidance/standards. In accordance with EN ISO 17664, these instructions provide two validated methods for cleaning: a manual and an automatic method (using an automatic washing machine). In accordance with ISO 17665, the half cycle "overkill" approach was used for sterilization validation and demonstrates a sterility assurance level (SAL) of 10⁻⁶.

ABANZA has validated the reprocessing instructions as capable and effective for reuse the components of the WasherCap™ Surgical Instrument Set. These instructions are described in more detail, in a document called "Reprocessing Instructions" (CI-001) supplied with the Surgical Instrument Set. However, it remains the responsibility of the responsible user to ensure that the reprocessing of the set is carried out correctly and effectively, in compliance with legislation and ordinances and using equipment and material approved by the FDA, CE marked or authorized by the local government body responsible for them.

D. PRE-CLEANING

- Point of use reprocessing**
 - Remove gross soil (liquid, tissue, bone fragments, organic debris, etc.) from the instruments, sterilization tray and lid with a disposable cloth or allowing water through them as soon as reasonably possible after use.
 - Do not allow blood, body fluids or other organic debris to dry on the instruments set.
- Transport to reprocessing area**
 - Take the instruments, sterilization tray and lid to the reprocessing area as soon as possible (covered or closed) to avoid contamination risks and avoiding bumping the instruments against each other during transportation.
 - It is recommended to put the instruments inside the sterilization tray during transport, to prevent them from bumping against each other.
- Preparation for cleaning in the reprocessing area (for manual and automatic)**
 - Remove the instruments from the sterilization tray before performing the cleaning process.
 - Disassemble the instruments consisting of several parts in the reprocessing area for their proper cleaning. Refer to the "Reprocessing Instructions – Appendix II" (CI-001) supplied with the instrument set to identify the instruments that must be disassembled before cleaning.

E. MANUAL CLEANING

STEPS:

- Washing**
 - If it has not been done previously at the point of use of the instruments, remove excess dirt from the instruments with a disposable cloth.
 - Prepare a cleaning solution using an enzymatic or alkaline cleaning agent (pH equal to or less than 12) in accordance with the manufacturer's instructions for its use in terms of recommended types of materials, concentrations, temperatures, water quality, exposure times and cleaning method.
 - Completely immerse the instruments, sterilization tray and lid separately (each individually and not inside the sterilization tray).
- Brushing**
 - Remove any residue or traces of blood found on instruments, sterilization tray and lid using a soft smooth bristles brush.
- Rinsing 1**
 - Rinse the instruments, sterilization tray and lid in running water (Service water in accordance with AAMI TIR34: Hardness < 150 ppm, Conductivity < 500 µS/cm, Chloride < 250 ppm, Total organic carbon < 1 mg/mL) until all visible traces of detergent have been removed.
- Ultrasonic Bath**
 - Prepare a solution with a cleaning agent suitable for the ultrasonic bath, in accordance with the manufacturer's instructions for the use of the cleaner in terms of recommended types of materials, concentrations, temperatures, water quality, exposure times and cleaning method.
 - Immerse the instruments, sterilization tray and lid in an ultrasonic bath.
 - Cleaning for a minimum of 10 minutes at a frequency of 40-50 kHz its recommended, provided that manufacturer's instructions advise so.
- Rinsing 2 (after ultrasonic bath)**
 - Rinse the instruments, sterilization tray and lid in running water (Utility water in accordance with AAMI TIR34: Hardness < 150 ppm, Conductivity < 500 µS/cm, Chlorides < 250 ppm, Total organic carbon < 1 mg/mL) until all visible traces of detergent have been removed.
- Final rinsing**
 - Rinse the instruments, sterilization tray and lid in DEIONIZED water (Critical water in accordance with AAMI 5 TIR34: Hardness < 1 ppm, Conductivity < 10 µS/cm, Chloride < 1 ppm, Bacteria < 10 CFU/mL, Endotoxins < 10 CFU/mL).
- Drying**
 - Dry the instruments, sterilization tray and lid with a clean cloth made of fiber-free materials to ensure that the parts are completely dry. Filter pressurized air can be used to help with drying.

F. DISINFECTION (Optional)

Under certain risk classifications, the World Health Organization (WHO) or local regulatory authorities recommend special Creutzfeldt-Jakob Disease (CJD/TSE - Transmissible Spongiform Encephalopathy) inactivation procedures. ABANZA reusable instruments are not used in surgical procedures where they come into contact with CJD/TSE infectious tissue (according to WHO). Therefore, excessively aggressive decontamination procedures (using aggressive detergents such as sodium hydroxide or sodium hypochlorite, aldehydes, bromides, iodides) are not required, and are not recommended for normal reprocessing because material degradation may occur. Refer to the WHO and local regulations for more information.

G. AUTOMATIC CLEANING

Note: To carry out automatic cleaning, automatic washing machines FDA approved (FDA or CE marked according to ISO 15883), properly installed, qualified and subject to regular maintenance and testing should be used. The operating instructions of the automatic washing machine manufacturer should be followed. Use only cleaning agents recommended by the manufacturer specifically for their use in automatic washing machines.

STEPS:

- Loading into the automatic washing machine**
 - Introduce the instruments, sterilization tray and lid separately (each individually and not inside the sterilization tray) into the washing machine, disassembled (Refer to the "Reprocessing Instructions – Appendix II" (CI-001) supplied with the instrument set to identify the instruments that must be disassembled before cleaning) and in an orderly manner.
 - Do not stack the instruments, sterilization tray and lid on top of each other.
- Cleaning cycle**
 - Process instruments, sterilization tray and lid with a standard cleaning cycle for surgical instruments validated by the automatic washing machine manufacturer using a detergent suitable for automatic cleaning.
 - The parameters of the cycle validated by ABANZA are listed below:

PARAMETERS OF THE AUTOMATIC WASHING CYCLE			
Phase	Temperature	Duration	Concentration detergent
Prewashing I	10 °C (50 °F)	10 minutes	NA
Washing	55 °C (131°F)	5 minutes	0.3 % detergent
Neutralization	10 °C (50 °F)	2 minutes	NA
Rinsing II	10 °C (50 °F)	1 minutes	NA

III. Remove the elements from the automatic washing machine

- Remove the instruments, sterilization tray and lid from the automatic washing machine and check if there is any remaining dirt; if so, repeat the process.

IV. Drying

- Dry the instruments, sterilization tray and lid with a clean cloth made of a fiber-free material to ensure that they are completely dry. Filter pressurized air can be used to help with drying.

H. INSPECTION

After cleaning process, the instruments, sterilization tray and lid should be inspected before preparing them for sterilization. For this purpose, each instrument, sterilization tray and lid should be visually inspected and checked for signs of visible dirt and/or damage and wear.

- In terms of cleaning, complete removal of dirt from all instruments, sterilization tray and lid surfaces should be ensured.
- In terms of damage and wear, it is considered unacceptable and therefore as failure criteria when any damage to the instruments, sterilization tray and lid is identified visually, including but not limited to: burrs, breakage, corrosion (rust, pitting), discoloration, excessive scratches, flaking, wear, cracks, malfunction, etc. should be visually checked. Also check that the references marked on the instruments, sterilization tray and lid are legible in order to maintain their traceability.

In the event of damage or wear that could affect the proper functioning of the instruments, sterilization tray and lid ask for a replacement to the ABANZA representative.

I. PREPARATION FOR STERILIZATION

STEPS:

- Reassembly the instruments that require disassembly during the cleaning process before the sterilization process in the reprocessing area. Refer to the "Reprocessing Instructions – Appendix II" (CI-001) supplied with the instrument set to identify the instruments that must be reassembled before sterilization.
- Place the instruments in the sterilization tray for sterilization.
- Close the tray with its lid and check that the latches are securely fastened to prevent them from opening.
- Place the closed sterilization tray in an authorized sterilization container (FDA approved or CE marked) with a tight lid and suitable for steam sterilization at the user's discretion.

J. STERILIZATION PROCESS

Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions. Refer to "Reprocessing Instructions" (CI-001) supplied with the instrument set for reprocessing.

- Perform sterilization of reusable surgical instruments by moist heat method in an autoclave (sterilizer) with a pre-vacuum cycle (forced air removal).

Note: The autoclave must comply with the requirements of EN 285, EN 13060, ISO 17665 and AINSI/AAMI ST79, and be validated and maintained accordingly. It is recommended to follow the parameters validated by ABANZA and described below.

MOIST HEAT STERILIZATION PARAMETERS:

STERILIZATION PARAMETERS: FOR THE USA ONLY				
Method	Sterilizer cycle	Temperature of exposure	Exposure time	Minimum drying time
Moist heat sterilization	Pre-vacuum	132 °C (270 °F)	4 minutes (minimum)	30 minutes (minimum)

STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY				
Method	Sterilizer cycle	Temperature of exposure	Exposure time	Minimum drying time
Moist heat sterilization	Pre-vacuum	132 °C to 137 °C (270 °F to 277 °F)	4 minutes (minimum)	30 minutes (minimum)

K. PRE-USE STORAGE

- Place the sterilization container with the sterile instruments in a specific area with limited access that is well ventilated and offers protection against dust, humidity, insects, temperature and humidity extremes.

Note: The set should be storage inside the sterilization container. If there is any indication or sign that the sterilization container has been opened or compromised, the instruments should be re-sterilized. Storage period depends on several factors (storage form, environmental conditions, handling, etc.), so the maximum storage period for sterilized reusable instruments before use should be defined by each user.

L. SPECIAL CAUTION: TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy Agents.

The agents for transmission of Creutzfeldt-Jakob disease are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing method of decontamination and sterilization as described above may not be appropriate where Creutzfeldt-Jakob transmission is a risk.

In general, the tissues that come into contact with orthopedic surgical instruments are those of low Transmissible Spongiform Encephalopathy (TSE) infectivity. However, particular precautions should be taken when handling instruments that have been used on known, suspected, or at risk-patients.

M. PRECAUTIONS

- The WasherCap™ Surgical Instrument Set should only be accepted if the original manufacturer's packaging and labelling are present and arrive intact. If the package has been opened and/or if any component of the Set is damaged, inform your local distributor or ABANZA representative for replacement.
- To avoid damaging the instruments that are designed to be turned or screwed must not be hit or blunted. When two devices are intended to be threaded together, ensure that they are fully engaged prior to use.
- Do not use ABANZA instruments for any purpose other than their intended use.
- Do not use an instrument that is intended to be used with a specific implant on another implant.
- Some instruments in the set, due to their intended use, have sharp edges. To avoid possible damage to the tendons, patient, and/or users, do not use the instruments for other than the indicated use.
- For the correct use of the **countersink** instrument, it is required:
 - Disassemble (main body and top element) before carrying out cleaning procedures at the reprocessing area, with the help of the countersink tool.
 - Reassemble it for sterilization at the reprocessing area. When assembling, it is important not to mix the parts between the different sizes of countersink (observe the size engraved on each instrument). Check that it has been correctly clamped (when fully threaded) and, if necessary, tighten it with the countersink tool.
 - For further information about the assembly and disassembly of the instruments, please refer to the "Reprocessing Instructions – Appendix II" (CI-001) provided by ABANZA with the Instrument Set.

N. WARNINGS

- The instruments of the WasherCap™ Surgical Instrument Set should be only handled by surgeons who are familiar with the clinical use, surgical technique, proper selection and placement of the WasherCap™ device and the use of it associated Surgical Instrument Set.
- Federal (US) law restricts this device to sale by or on the order of a physician.
- The instruments of the WasherCap™ Surgical Instrument Set are supplied non-sterile and must be reprocessed (cleaned and sterilized) before use. For this purpose, it is recommended to follow the instructions cited in these Instructions for Use and in the document "Reprocessing Instructions" (CI-001) provided by ABANZA with the Surgical Instrument Set.
- Do not use the instruments for other purposes or surgical procedures other than those indicated in these Instructions for Use.

O. DISCLAIMER

Please consult the Instructions For Use (IFU) of the associated device to this Surgical Instrument Set for further information. Please contact your local distributor or ABANZA representative for a copy of the IFU of the associated device or the "Reprocessing Instructions" (CI-001) of the WasherCap™ Surgical Instrument Set.

All surgical instruments which may be damaged, scratched, showing defects, corrosion, cracked joints or discoloration, should not be used and be managed according to hospital protocol rules. Any misuse, modification or replacement of the instruments is forbidden. If this recommendation is not followed, ABANZA rejects any responsibility derived from any subsequent consequence.

This document is not subject to control after distribution, and thus, it is the responsibility of the reader to confirm that this is the most current revision or to contact ABANZA TECNOMED's customer service to request the latest version.

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